



ETO Sterilizer:

**Sterilize your equipments with our top
quality products**

- ❖ **ETO Gas Sterilizer**
- ❖ **ETO Sterilization Plants**
- ❖ **High Pressure Steam Sterilizer**
- ❖ **Pharmaceutical ETO Sterilizers**
- ❖ **Pharma Autoclave - Medical Autoclave**



ETO Gas Sterilizer

The ETO (Ethylene Oxide) sterilization is the process of sterilizing the products that are used in medical and pharmaceuticals. The ETO gas sterilization process is carried out in four different constituents:

the concentration of gas, humidity, temperature and time required to carry out the process of sterilization. We offer a wide assortment of ETO gas sterilizer that work effectively and can handle high temperature and pressure.





ETO Sterilization Plants

The sterilization process of all the products is carried out in special plants. These ETO sterilization plants are made to execute 3 main procedures: pre-conditioning, sterilization and aeration. This process needs proper isolation to make completely sterilized products. They are offered in various size and specifications.





High Pressure Steam Sterilizer

The sterilization process needs high temperature and pressure to effectively sterilize the medical and pharmaceutical products. We offer high pressure steam sterilizer that are used to sanitize the product so that they are completely safe for use in medical purpose. They are effective in cleaning the hospital appliances from microorganisms.





Pharmaceutical ETO Sterilizers

The ETO sterilizers are used in various industries to sanitize their products and make them completely safe for use. We are the top manufacturer of pharmaceutical ETO sterilizers that are highly demanded in pharmaceutical industries for their tools, machines and devices to make hygienic and safe medicines.





Pharma - Medical Autoclave

When the steam sterilization process is carried out in weight chamber that is called an autoclave. We offer wide range of pharma autoclave and medical autoclave that are used for cleaning infusion needles, utensils and surgical devices that are used in medicals. We offer horizontal high pressure autoclave that can execute the process at high temperature and pressure.





Contact Us

Address: No. 7, Khodiyar Estate, Near Sakariba
Estate, Behind Indo German Tool Room,
Phase 4, G. I. D. C., Vatva,
Ahmedabad - 382445, Gujarat, India

Email: info@krishna-engg.com